

Product Datasheet


FECH Antibody (orb1241847)

Description

FECH Antibody


Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human FECH.
Target	FECH
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store FECH antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	FECH antibody can be used for detection of FECH by ELISA at 1:1562500. FECH antibody can be used for detection of FECH by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	47 kDa, 48 kDa
Uniprot ID	Q8NANO
NCBI	NP_001012533
Dilution Range	FECH antibody can be used for detection of FECH by ELISA at 1:1562500. FECH antibody can be used for detection of FECH by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

FECH



See Immunoblot Data and Customer Feedback for more information

Antibody used in WB on Mouse liver, hear...



Antibody used in WB on Human Jurkat cell...