

Product Datasheet

CYP2D6 Antibody (orb1241849)

Description

CYP2D6 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Antibody
used in WB
on Human
Liver at 2....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CYP2D6.

Target

CYP2D6

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance


Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CYP2D6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.



Antibody
used in IHC
on Human
Kidney at ...

Note

For research use only

Application notes

CYP2D6 antibody can be used for detection of CYP2D6 by ELISA at 1:312500. CYP2D6 antibody can be used for detection of CYP2D6 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

55 kDa

Uniprot ID

P10635

NCBI

NP_000097

Dilution Range

CYP2D6 antibody can be used for detection of CYP2D6 by ELISA at 1:312500. CYP2D6 antibody can be used for detection of CYP2D6 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.