

Product Datasheet

GGTLC1 Antibody (orb1241858)

Description

GGTLC1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GGTLA4.

Target

GGTLC1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GGTLA4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

GGTLA4 antibody can be used for detection of GGTLA4 by ELISA at 1:312500. GGTLA4 antibody can be used for detection of GGTLA4 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

25 kDa, 61 kDa, 61 kDa, 61 kDa, 61 kDa, 24 kDa, 24 kDa, 27 kDa, 24 kDa, 27 kDa, 24 kDa, 62 kDa, 62 kDa, 62 kDa

Uniprot ID

[Q9BX51](#)

NCBI

[NP_563577](#)

Dilution Range

GGTLA4 antibody can be used for detection of GGTLA4 by ELISA at 1:312500. GGTLA4 antibody can be used for detection of GGTLA4 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.


90kDa

60kDa

30kDa

20kDa

Antibody used in WB on Human kidney at 5...



Antibody used in IHC on Human Kidney at ...