

Product Datasheet

SLC22A1 Antibody (orb1241873)

Description

SLC22A1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC22A1.

Target

SLC22A1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC22A1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC22A1 antibody can be used for detection of SLC22A1 by ELISA at 1:62500. SLC22A1 antibody can be used for detection of SLC22A1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

61 kDa, 56 kDa

Uniprot ID

[O15245](#)

NCBI

[NP_003048](#)

Dilution Range


SLC22A1 antibody can be used for detection of SLC22A1 by ELISA at 1:62500. SLC22A1 antibody can be used for detection of SLC22A1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.

70 kDa
60 kDa
50 kDa
40 kDa
30 kDa
20 kDa

Antibody
used in WB
on Human
Jurkat 2.5 ...



Antibody
used in IHC
on Human
Intestine ...



Antibody
used in WB
on Jurkat 2.5
ug/ml ...