

Product Datasheet

SLC13A3 Antibody (orb1241890)


Description

SLC13A3 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC13A3.
Target	SLC13A3
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC13A3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SLC13A3 antibody can be used for detection of SLC13A3 by ELISA at 1:312500. SLC13A3 antibody can be used for detection of SLC13A3 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	61 kDa, 67 kDa
Uniprot ID	Q8WWT9
NCBI	NP_001011554
Dilution Range	SLC13A3 antibody can be used for detection of SLC13A3 by ELISA at 1:312500. SLC13A3 antibody can be used for detection of SLC13A3 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

80L
40L
20L

Antibody used in WB on Human Jurkat 2.5 ...



Antibody used in IHC on Human Kidney at ...