

Product Datasheet

MAT1A Antibody (orb1241903)


Description

MAT1A Antibody

Species/Host	Rabbit
Reactivity	C. elegans, Canine, Drosophila, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MAT1A.
Target	MAT1A
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store MAT1A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	MAT1A antibody can be used for detection of MAT1A by ELISA at 1:1562500. MAT1A antibody can be used for detection of MAT1A by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	Q00266
NCBI	NP_000420
Dilution Range	MAT1A antibody can be used for detection of MAT1A by ELISA at 1:1562500. MAT1A antibody can be used for detection of MAT1A by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Jurkat 1.25...



Antibody used in IHC on Human Kidney at ...