

Product Datasheet

SLC12A1 Antibody (orb1241905)

Description

SLC12A1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC12A1.

Target

SLC12A1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC12A1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC12A1 antibody can be used for detection of SLC12A1 by ELISA at 1:312500. SLC12A1 antibody can be used for detection of SLC12A1 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

47 kDa, 121 kDa

Uniprot ID

[Q8IUN5](#)

NCBI

[AAH40138](#)

Dilution Range


SLC12A1 antibody can be used for detection of SLC12A1 by ELISA at 1:312500. SLC12A1 antibody can be used for detection of SLC12A1 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human Jurkat 0.2-...



Antibody used in IHC on Human Kidney at ...