


Product Datasheet

GNB1L Antibody (orb1241911)


Description

GNB1L Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GNB1L.
Target	GNB1L
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GNB1L antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	GNB1L antibody can be used for detection of GNB1L by ELISA at 1:12500. GNB1L antibody can be used for detection of GNB1L by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	Q9BYB4
NCBI	NP_443730
Dilution Range	GNB1L antibody can be used for detection of GNB1L by ELISA at 1:12500. GNB1L antibody can be used for detection of GNB1L by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Skin at 4....



Antibody used in WB on Human Jurkat 2.5 ...



Antibody used in WB on Jurkat 2.5 ug/ml ...