


Product Datasheet

WNT2B Antibody (orb1241922)


Description

WNT2B Antibody


Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human WNT2B.
Target	WNT2B
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store WNT2B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	WNT2B antibody can be used for detection of WNT2B by ELISA at 1:62500. WNT2B antibody can be used for detection of WNT2B by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	41 kDa, 44 kDa
Uniprot ID	Q5T6U1
NCBI	NP_004176
Dilution Range	WNT2B antibody can be used for detection of WNT2B by ELISA at 1:62500. WNT2B antibody can be used for detection of WNT2B by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Lung at 4....



Antibody used in IHC on Human Testis.



Antibody used in WB on Jurkat 1.25 ug/ml...