



Product Datasheet

EXOSC6 Antibody (orb1241975)

Description	EXOSC6 Antibody	
Species/Host	Rabbit	
Reactivity	Canine, Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human EXOSC6.	Antibody used in IHC on Human Muscle at ...
Target	EXOSC6	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store EXOSC6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in IHC on Human Liver at 4...
Note	For research use only	
Application notes	EXOSC6 antibody can be used for detection of EXOSC6 by ELISA at 1:312500. EXOSC6 antibody can be used for detection of EXOSC6 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	30 kDa	
Uniprot ID	Q5RKV6	Antibody used in WB on Human Jurkat 5.0 ...
NCBI	NP_478126	
Dilution Range	EXOSC6 antibody can be used for detection of EXOSC6 by ELISA at 1:312500. EXOSC6 antibody can be used for detection of EXOSC6 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	