

Product Datasheet

NXF5 Antibody (orb1241981)

Description

NXF5 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NXF5.

Target

NXF5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store NXF5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

NXF5 antibody can be used for detection of NXF5 by ELISA at 1:312500. NXF5 antibody can be used for detection of NXF5 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

44 kDa, 72 kDa, 72 kDa, 72 kDa, 72 kDa, 72 kDa

Uniprot ID

[Q9H1B4](#)

NCBI


[NP_116564](#)

Dilution Range


NXF5 antibody can be used for detection of NXF5 by ELISA at 1:312500. NXF5 antibody can be used for detection of NXF5 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Kidney at ...



Antibody used in IHC on Human Intestine ...



Antibody used in WB on Human Jurkat 2.5 ...