

Product Datasheet

NOC4L Antibody (orb1242001)


Description

NOC4L Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NOC4L.
Target	NOC4L
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store NOC4L antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	NOC4L antibody can be used for detection of NOC4L by ELISA at 1:12500. NOC4L antibody can be used for detection of NOC4L by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	57 kDa, 58 kDa
Uniprot ID	Q9BVI4
NCBI	NP_076983
Dilution Range	NOC4L antibody can be used for detection of NOC4L by ELISA at 1:12500. NOC4L antibody can be used for detection of NOC4L by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Jurkat 0.2-...



Antibody used in IHC on HUMAN Kidney at ...