

Product Datasheet

DAZAP1 Antibody (orb1242014)

Description

DAZAP1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human DAZAP1.

Target

DAZAP1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store DAZAP1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

DAZAP1 antibody can be used for detection of DAZAP1 by ELISA at 1:312500. DAZAP1 antibody can be used for detection of DAZAP1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

45 kDa

Uniprot ID

[Q96EP5](#)

NCBI


[NP_061832](#)

Dilution Range


DAZAP1 antibody can be used for detection of DAZAP1 by ELISA at 1:312500. DAZAP1 antibody can be used for detection of DAZAP1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Liver at 4...



Antibody used in IHC on Human Heart at 4...



Antibody used in WB on Human Daudi at 1....