


Product Datasheet

MATR3 Antibody (orb1242017)


Description

MATR3 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MATR3.
Target	MATR3
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store MATR3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	MATR3 antibody can be used for detection of MATR3 by ELISA at 1:312500. MATR3 antibody can be used for detection of MATR3 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	93 kDa, 95 kDa
Uniprot ID	P43243
NCBI	NP_061322
Dilution Range	MATR3 antibody can be used for detection of MATR3 by ELISA at 1:312500. MATR3 antibody can be used for detection of MATR3 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Heart at 4...



Antibody used in IHC on Human Muscle at ...



Antibody used in WB on Human Jurkat cell...