

Product Datasheet

LARP7 Antibody (orb1242030)

Description

LARP7 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Jurkat 1.25...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human LARP7.

Target

LARP7

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store LARP7 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

LARP7 antibody can be used for detection of LARP7 by ELISA at 1:312500. LARP7 antibody can be used for detection of LARP7 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

39 kDa, 67 kDa, 67 kDa

Uniprot ID

[Q4G0J3](#)

NCBI

[EAX06285](#)

Dilution Range

LARP7 antibody can be used for detection of LARP7 by ELISA at 1:312500. LARP7 antibody can be used for detection of LARP7 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.