

Product Datasheet

TRNT1 Antibody (orb1242038)

Description

TRNT1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Antibody
used in WB
on Human
brain at 2....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TRNT1.

Target

TRNT1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance


Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TRNT1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.



Antibody
used in IHC
on Human
Kidney at ...

Note

For research use only

Application notes

TRNT1 antibody can be used for detection of TRNT1 by ELISA at 1:62500. TRNT1 antibody can be used for detection of TRNT1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

45 kDa, 48 kDa

Uniprot ID

[Q96Q11](#)

NCBI

[NP_057084](#)

Dilution Range

TRNT1 antibody can be used for detection of TRNT1 by ELISA at 1:62500. TRNT1 antibody can be used for detection of TRNT1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.