

Product Datasheet

CSTF2T Antibody (orb1242041)

Description

CSTF2T Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Brain at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CSTF2T.

Target

CSTF2T

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CSTF2T antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CSTF2T antibody can be used for detection of CSTF2T by ELISA at 1:1562500. CSTF2T antibody can be used for detection of CSTF2T by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

64 kDa

Uniprot ID

[Q9H0L4](#)

NCBI

[NP_056050](#)

Dilution Range

CSTF2T antibody can be used for detection of CSTF2T by ELISA at 1:1562500. CSTF2T antibody can be used for detection of CSTF2T by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.