

Product Datasheet

SNRNP35 Antibody (orb1242064)

Description

SNRNP35 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human U1SNRNPBP.

Antibody used in WB on Human Jurkat 2.5 ...

Target

SNRNP35

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store U1SNRNPBP antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in IHC on Human Kidney at ...

Note

For research use only

Application notes

U1SNRNPBP antibody can be used for detection of U1SNRNPBP by ELISA at 1:1562500. U1SNRNPBP antibody can be used for detection of U1SNRNPBP by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

28 kDa, 29 kDa, 29 kDa, 30 kDa

Uniprot ID

[Q5XKN9](#)

NCBI

[NP_851034](#)

Dilution Range

U1SNRNPBP antibody can be used for detection of U1SNRNPBP by ELISA at 1:1562500. U1SNRNPBP antibody can be used for detection of U1SNRNPBP by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

80kDa

60kDa

50kDa

40kDa

