

Product Datasheet

NUDT21 Antibody (orb1242066)


Description

NUDT21 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NUDT21.
Target	NUDT21
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store NUDT21 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	NUDT21 antibody can be used for detection of NUDT21 by ELISA at 1:312500. NUDT21 antibody can be used for detection of NUDT21 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	25 kDa
Uniprot ID	O43809
NCBI	NP_008937
Dilution Range	NUDT21 antibody can be used for detection of NUDT21 by ELISA at 1:312500. NUDT21 antibody can be used for detection of NUDT21 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HepG2 at 1....



Antibody used in IHC on Human Intestine ...