

Product Datasheet

RNASEH2A Antibody (orb1242086)

Description


RNASEH2A Antibody

Species/Host	Rabbit
Reactivity	Human, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RNASEH2A.
Target	RNASEH2A
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RNASEH2A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	RNASEH2A antibody can be used for detection of RNASEH2A by ELISA at 1:312500. RNASEH2A antibody can be used for detection of RNASEH2A by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	O75792
NCBI	NP_006388
Dilution Range	RNASEH2A antibody can be used for detection of RNASEH2A by ELISA at 1:312500. RNASEH2A antibody can be used for detection of RNASEH2A by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.






Antibody used in WB on Human Jurkat 1.25...



Antibody used in IHC on Human Kidney at ...