

Product Datasheet

NOP56 Antibody (orb1242089)

Description

NOP56 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NOL5A.

Target

NOP56

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store NOL5A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

NOL5A antibody can be used for detection of NOL5A by ELISA at 1:62500. NOL5A antibody can be used for detection of NOL5A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

66 kDa

Uniprot ID

[O00567](#)

NCBI

[NP_006383](#)

Dilution Range

NOL5A antibody can be used for detection of NOL5A by ELISA at 1:62500. NOL5A antibody can be used for detection of NOL5A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.


90 kDa

40 kDa

31 kDa

22 kDa

Antibody used in WB on Human 721_B at 1 ...



Antibody used in IHC on Human Placenta a...