



Product Datasheet

NXF1 Antibody (orb1242092)


Description

NXF1 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NXF1.
Target	NXF1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store NXF1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	NXF1 antibody can be used for detection of NXF1 by ELISA at 1:312500. NXF1 antibody can be used for detection of NXF1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	68 kDa, 40 kDa
Uniprot ID	Q9UBU9
NCBI	NP_006353
Dilution Range	NXF1 antibody can be used for detection of NXF1 by ELISA at 1:312500. NXF1 antibody can be used for detection of NXF1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human kidney at ...



Antibody used in WB on Human HepG2 at 1....



See Immunoblot 2 Data and Customer Feedback for more information

Antibody used in WB on Human HeLa at: 1:...