



Product Datasheet

NFKB1 antibody (orb12421)


Description	Rabbit polyclonal antibody to NFKB1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	/
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	120KDa
Uniprot ID	P19838
Entrez	4790
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200, IF/ICC: 1:20-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF7 cell...



Immunohistochemical staining of human ga...



Immunohistochemical staining of human to...