

Product Datasheet

GNB2L1 Antibody (orb1242102)

Description

GNB2L1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GNB2L1.

Target

GNB2L1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GNB2L1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

GNB2L1 antibody can be used for detection of GNB2L1 by ELISA at 1:62500. GNB2L1 antibody can be used for detection of GNB2L1 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

35 kDa

Uniprot ID

[P63244](#)

NCBI

[NP_006089](#)

Dilution Range

GNB2L1 antibody can be used for detection of GNB2L1 by ELISA at 1:62500. GNB2L1 antibody can be used for detection of GNB2L1 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

GNB2L1

40 kDa
30 kDa
20 kDa

Antibody used in WB on Human Jurkat 0.2-...

GNB2L1

40 kDa
30 kDa
20 kDa

Antibody used in WB on Human MCF-7 at 1....