


Product Datasheet

THOC1 Antibody (orb1242115)


Description

THOC1 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human THOC1.
Target	THOC1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store THOC1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	THOC1 antibody can be used for detection of THOC1 by ELISA at 1:12500. THOC1 antibody can be used for detection of THOC1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	76 kDa
Uniprot ID	Q96FV9
NCBI	NP_005122
Dilution Range	THOC1 antibody can be used for detection of THOC1 by ELISA at 1:12500. THOC1 antibody can be used for detection of THOC1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Kidney at ...



Antibody used in IHC on Human Lung at 4....



Antibody used in IHC on human prostate t...