

Product Datasheet

ZRSR2 Antibody (orb1242117)


Description

ZRSR2 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZRSR2.
Target	ZRSR2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ZRSR2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	ZRSR2 antibody can be used for detection of ZRSR2 by ELISA at 1:62500. ZRSR2 antibody can be used for detection of ZRSR2 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	53 kDa
Uniprot ID	Q15696
NCBI	NP_005080
Dilution Range	ZRSR2 antibody can be used for detection of ZRSR2 by ELISA at 1:62500. ZRSR2 antibody can be used for detection of ZRSR2 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HepG2 at 2....



Antibody used in IHC on Human Kidney at ...