


Product Datasheet

SFPQ Antibody (orb1242118)


Description

SFPQ Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SFPQ.
Target	SFPQ
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SFPQ antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SFPQ antibody can be used for detection of SFPQ by ELISA at 1:12500. SFPQ antibody can be used for detection of SFPQ by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	78 kDa
Uniprot ID	P23246
NCBI	NP_005057
Dilution Range	SFPQ antibody can be used for detection of SFPQ by ELISA at 1:12500. SFPQ antibody can be used for detection of SFPQ by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Intestine ...



Antibody used in WB on Human HepG2 at 1....



Antibody used in IHC on Human Skin at 4....