

Product Datasheet

STAU1 Antibody (orb1242127)


Description

STAU1 Antibody


Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human STAU1.
Target	STAU1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store STAU1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	STAU1 antibody can be used for detection of STAU1 by ELISA at 1:62500. STAU1 antibody can be used for detection of STAU1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	63 kDa, 56 kDa, 55 kDa, 55 kDa, 55 kDa
Uniprot ID	O95793
NCBI	NP_059347
Dilution Range	STAU1 antibody can be used for detection of STAU1 by ELISA at 1:62500. STAU1 antibody can be used for detection of STAU1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HepG2 at 1....



Antibody used in IHC on Human Intestine ...



Antibody used in WB on Human HepG2 cells...