

Product Datasheet

PABPC4 Antibody (orb1242142)


Description

PABPC4 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PABPC4.
Target	PABPC4
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PABPC4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PABPC4 antibody can be used for detection of PABPC4 by ELISA at 1:312500. PABPC4 antibody can be used for detection of PABPC4 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	71 kDa
Uniprot ID	Q13310
NCBI	NP_003810
Dilution Range	PABPC4 antibody can be used for detection of PABPC4 by ELISA at 1:312500. PABPC4 antibody can be used for detection of PABPC4 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HepG2 at 1....



Antibody used in IHC on Human Kidney at ...