

Product Datasheet

SNRPF Antibody (orb1242160)

Description


SNRPF Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SNRPF.
Target	SNRPF
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SNRPF antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SNRPF antibody can be used for detection of SNRPF by ELISA at 1:312500. SNRPF antibody can be used for detection of SNRPF by western blot at 0.31 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	9 kDa
Uniprot ID	P62306
NCBI	NP_003086
Dilution Range	SNRPF antibody can be used for detection of SNRPF by ELISA at 1:312500. SNRPF antibody can be used for detection of SNRPF by western blot at 0.31 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HepG2 at 1

...



Antibody used in IHC on Human Heart at 4...