

Product Datasheet

RBMS1 Antibody (orb1242165)

Description


RBMS1 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RBMS1.
Target	RBMS1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RBMS1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	RBMS1 antibody can be used for detection of RBMS1 by ELISA at 1:1562500. RBMS1 antibody can be used for detection of RBMS1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	44 kDa, 45 kDa, 42 kDa
Uniprot ID	P29558
NCBI	NP_002888
Dilution Range	RBMS1 antibody can be used for detection of RBMS1 by ELISA at 1:1562500. RBMS1 antibody can be used for detection of RBMS1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.






Antibody used in WB on Human Jurkat 2.5 ...



Antibody used in IHC on Human Skin at 4....