

Product Datasheet


RBM4 Antibody (orb1242166)

Description


RBM4 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RBM4.
Target	RBM4
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RBM4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	RBM4 antibody can be used for detection of RBM4 by ELISA at 1:62500. RBM4 antibody can be used for detection of RBM4 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	Q96LT9
NCBI	NP_002887
Dilution Range	RBM4 antibody can be used for detection of RBM4 by ELISA at 1:62500. RBM4 antibody can be used for detection of RBM4 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

RBM4



Antibody used in WB on Hum. Adult Placen...



Antibody used in WB on Human MCF7 cells ...