

Product Datasheet

SF3B1 Antibody (orb1242230)


Description

SF3B1 Antibody


Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SF3B1.
Target	SF3B1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SF3B1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SF3B1 antibody can be used for detection of SF3B1 by ELISA at 1:62500. SF3B1 antibody can be used for detection of SF3B1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	143 kDa, 16 kDa
Uniprot ID	O75533
NCBI	NP_036565
Dilution Range	SF3B1 antibody can be used for detection of SF3B1 by ELISA at 1:62500. SF3B1 antibody can be used for detection of SF3B1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Thymus at 1...



Antibody used in IHC on Human Liver at 4...



Antibody used in IHC on Human Skin at 16...