

Product Datasheet

DAZ4 Antibody (orb1242232)

Description

DAZ4 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human DAZ4.

Target

DAZ4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store DAZ4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

DAZ4 antibody can be used for detection of DAZ4 by ELISA at 1:312500. DAZ4 antibody can be used for detection of DAZ4 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

36 kDa, 65 kDa, 60 kDa, 41 kDa, 33 kDa, 65 kDa, 63 kDa, 50 kDa, 44 kDa

Uniprot ID

[Q86SG3](#)

NCBI


[AAF91331](#)

Dilution Range


DAZ4 antibody can be used for detection of DAZ4 by ELISA at 1:312500. DAZ4 antibody can be used for detection of DAZ4 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Kidney at ...



Antibody used in IHC on Human Stomach at...



Antibody used in WB on Human HepG2 at 0....