

Product Datasheet

EXOSC3 Antibody (orb1242238)

Description

EXOSC3 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human EXOSC3.

Target

EXOSC3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store EXOSC3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

EXOSC3 antibody can be used for detection of EXOSC3 by ELISA at 1:312500. EXOSC3 antibody can be used for detection of EXOSC3 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

30 kDa, 17 kDa

Uniprot ID

[Q9NQT5](#)

NCBI


[NP_057126](#)

Dilution Range


EXOSC3 antibody can be used for detection of EXOSC3 by ELISA at 1:312500. EXOSC3 antibody can be used for detection of EXOSC3 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human Jurkat 0.2-...



Antibody used in IHC on Human Intestine ...