

## Product Datasheet

### ANP32A Antibody (orb1242254)

## Description

ANP32A Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, IHC, WB

### Applications

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ANP32A.

### Target

ANP32A

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ANP32A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

ANP32A antibody can be used for detection of ANP32A by ELISA at 1:312500. ANP32A antibody can be used for detection of ANP32A by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

29 kDa

### Uniprot ID

[P39687](#)

### NCBI

[NP\\_006296](#)


### Dilution Range

ANP32A antibody can be used for detection of ANP32A by ELISA at 1:312500. ANP32A antibody can be used for detection of ANP32A by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.

ANP32A




See Immunoblot 2 Data and Customer Feedback for more information

Antibody used in WB on Mouse cortex brai...



Antibody used in WB on Human Jurkat 0.2-...



Antibody used in IHC on human kidney at ...