

Product Datasheet

PRMT2 Antibody (orb1242257)


Description

PRMT2 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PRMT2.
Target	PRMT2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PRMT2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PRMT2 antibody can be used for detection of PRMT2 by ELISA at 1:1562500. PRMT2 antibody can be used for detection of PRMT2 by western blot at 0.625 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	48 kDa, 49 kDa
Uniprot ID	P55345
NCBI	NP_001526
Dilution Range	PRMT2 antibody can be used for detection of PRMT2 by ELISA at 1:1562500. PRMT2 antibody can be used for detection of PRMT2 by western blot at 0.625 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Jurkat 0.62...



Antibody used in WB on Human, Mouse at 1...



Antibody used in IHC on Human liver at 4...