

Product Datasheet

ZSCAN21 Antibody (orb1242313)

Description

ZSCAN21 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on
Transfected
293T ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZSCAN21.

Target

ZSCAN21

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZSCAN21 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZSCAN21 antibody can be used for detection of ZSCAN21 by ELISA at 1:62500. ZSCAN21 antibody can be used for detection of ZSCAN21 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

54 kDa

Uniprot ID

[Q9Y5A6](#)

NCBI

[NP_666019](#)

Dilution Range

ZSCAN21 antibody can be used for detection of ZSCAN21 by ELISA at 1:62500. ZSCAN21 antibody can be used for detection of ZSCAN21 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.