

Product Datasheet

ZNF296 Antibody (orb1242318)

Description

ZNF296 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF342.

Target

ZNF296

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF342 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZNF342 antibody can be used for detection of ZNF342 by ELISA at 1:62500. ZNF342 antibody can be used for detection of ZNF342 by western blot at 0.125 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

51 kDa

Uniprot ID

[Q8WUU4](#)

NCBI


[NP_660331](#)

Dilution Range


ZNF342 antibody can be used for detection of ZNF342 by ELISA at 1:62500. ZNF342 antibody can be used for detection of ZNF342 by western blot at 0.125 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Stomach.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human Jurkat 0.12...