

Product Datasheet

MINA Antibody (orb1242343)

Description

MINA Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MINA.

Target

MINA

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store MINA antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

MINA antibody can be used for detection of MINA by ELISA at 1:1562500. MINA antibody can be used for detection of MINA by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

32 kDa, 53 kDa, 53 kDa

Uniprot ID

[Q8IU6](#)

NCBI


[NP_694822](#)

Dilution Range

MINA antibody can be used for detection of MINA by ELISA at 1:1562500. MINA antibody can be used for detection of MINA by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Intestine.



Antibody used in WB on Human Thymus at 0...