

Product Datasheet

ZNF655 Antibody (orb1242370)

Description

ZNF655 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF655.

Target

ZNF655

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF655 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZNF655 antibody can be used for detection of ZNF655 by ELISA at 1:62500. ZNF655 antibody can be used for detection of ZNF655 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

54 kDa, 12 kDa, 20 kDa, 57 kDa

Uniprot ID

[Q8N720](#)

NCBI


[NP_001009960](#)

Dilution Range

ZNF655 antibody can be used for detection of ZNF655 by ELISA at 1:62500. ZNF655 antibody can be used for detection of ZNF655 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human HepG2 at 1....