

Product Datasheet

TEAD1 Antibody (orb1242379)

Description

TEAD1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human TEAD1.

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human TEAD1.

Target

TEAD1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TEAD1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TEAD1 antibody can be used for detection of TEAD1 by ELISA at 1:12500. TEAD1 antibody can be used for detection of TEAD1 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

46 kDa

Uniprot ID

[P28347](#)

NCBI

[NP_068780](#)

Dilution Range

TEAD1 antibody can be used for detection of TEAD1 by ELISA at 1:12500. TEAD1 antibody can be used for detection of TEAD1 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

95 kDa

46 kDa

25 kDa

22 kDa

Antibody used in WB on Human Placenta at...

TEAD1

90

65

40

25

22

Antibody Concentration: 1 µg/mL

Antibody used in WB on Human 721_B at 1 ...