

## Product Datasheet

### FEZF2 Antibody (orb1242412)

## Description

## FEZF2 Antibody

### Species/Host

Rabbit

### Reactivity

Canine, Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF312.

### Target

FEZF2

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF312 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

ZNF312 antibody can be used for detection of ZNF312 by ELISA at 1:1562500. ZNF312 antibody can be used for detection of ZNF312 by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

50 kDa, 52 kDa

### Uniprot ID

[Q8TBJ5](#)

### NCBI


[NP\\_060478](#)

### Dilution Range

ZNF312 antibody can be used for detection of ZNF312 by ELISA at 1:1562500. ZNF312 antibody can be used for detection of ZNF312 by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.



Antibody used in WB on Human HepG2 at 5....