

Product Datasheet

ZFR Antibody (orb1242426)

Description

ZFR Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Brain at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZFR.

Target

ZFR

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZFR antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZFR antibody can be used for detection of ZFR by ELISA at 1:312500. ZFR antibody can be used for detection of ZFR by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

67 kDa, 117 kDa

Uniprot ID

[Q6P5A3](#)

NCBI

[AAH62986](#)

Dilution Range

ZFR antibody can be used for detection of ZFR by ELISA at 1:312500. ZFR antibody can be used for detection of ZFR by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.