

Product Datasheet

PPARG Antibody (orb1242428)

Description

PPARG Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PPARG.

Target

PPARG

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PPARG antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PPARG antibody can be used for detection of PPARG by ELISA at 1:1562500. PPARG antibody can be used for detection of PPARG by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

56 kDa

Uniprot ID

[P37231](#)

NCBI


[NP_056953](#)

Dilution Range

PPARG antibody can be used for detection of PPARG by ELISA at 1:1562500. PPARG antibody can be used for detection of PPARG by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human Jurkat 1.25...