

Product Datasheet

RYBP Antibody (orb1242467)

Description

RYBP Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RYBP.

Target

RYBP

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store RYBP antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

RYBP antibody can be used for detection of RYBP by ELISA at 1:62500. RYBP antibody can be used for detection of RYBP by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

25 kDa

Uniprot ID

[Q8N488](#)

NCBI

[NP_036366](#)

Dilution Range

RYBP antibody can be used for detection of RYBP by ELISA at 1:62500. RYBP antibody can be used for detection of RYBP by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90kDa


40kDa

30kDa

20kDa

Antibody used in WB on Human Muscle at 0...

RYBP



See Remanded 2 Data and Customer Feedback for more information

Antibody used in WB on Human, Mouse at 1...