


Product Datasheet

ILF3 Antibody (orb1242469)

Description


ILF3 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ILF3.
Target	ILF3
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ILF3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	ILF3 antibody can be used for detection of ILF3 by ELISA at 1:62500. ILF3 antibody can be used for detection of ILF3 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	76 kDa, 95 kDa, 75 kDa
Uniprot ID	Q12906
NCBI	NP_004507
Dilution Range	ILF3 antibody can be used for detection of ILF3 by ELISA at 1:62500. ILF3 antibody can be used for detection of ILF3 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.




Antibody used in IHC on Human Lung.

ILF3



See Immunoblotting and Western Blotting for more information

Antibody used in WB on rat aorta at 1 ug...



Antibody used in WB on Human K562 cells ...