

Product Datasheet

DUX5 Antibody (orb1242473)

Description

DUX5 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human DUX5.
Target	DUX5
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store DUX5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	DUX5 antibody can be used for detection of DUX5 by ELISA at 1:12500. DUX5 antibody can be used for detection of DUX5 by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	19 kDa, 19 kDa
Uniprot ID	O43812
NCBI	NP_036281
Dilution Range	DUX5 antibody can be used for detection of DUX5 by ELISA at 1:12500. DUX5 antibody can be used for detection of DUX5 by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HepG2 at 5....