


Product Datasheet

SOX15 Antibody (orb1242491)

Description

SOX15 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SOX15.
Target	SOX15
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SOX15 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SOX15 antibody can be used for detection of SOX15 by ELISA at 1:312500. SOX15 antibody can be used for detection of SOX15 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	25 kDa
Uniprot ID	O60248
NCBI	NP_008873
Dilution Range	SOX15 antibody can be used for detection of SOX15 by ELISA at 1:312500. SOX15 antibody can be used for detection of SOX15 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in WB on Transfected 293T ...