

Product Datasheet

TADA3 Antibody (orb1242512)

Description

TADA3 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on
Transfected
293T ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TADA3L.

Target

TADA3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TADA3L antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TADA3L antibody can be used for detection of TADA3L by ELISA at 1:1562500. TADA3L antibody can be used for detection of TADA3L by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

49 kDa, 0 kDa

Uniprot ID

[O75528](#)

NCBI

[NP_006345](#)

Dilution Range

TADA3L antibody can be used for detection of TADA3L by ELISA at 1:1562500. TADA3L antibody can be used for detection of TADA3L by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

